Gebrauchsanweisung
Instruction manual
Manual de instrucciones

Rührwerk
Laboratory stirrer
Agitador de laboratorio
R100 C

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| | |
|---|--------------|
| Laboratory Stirrer R100 C Introduction | 34 |
| 1. General information | 35-36 |
| 1.1 Unpacking the instrument | 35 |
| 1.2 Scope of delivery | 35 |
| 1.3 Accessories (optional) | 36 |
| 1.4 Important instructions | 36 |
| 2. Safety Instructions | 37-38 |
| 2.1 Description of safety symbols in the instruction manual | 37 |
| 2.2 Description of safety symbols on the instrument | 37 |
| 2.3 General safety instructions | 37, 38 |
| 3. Mounting | 39-40 |
| 3.1 Mounting of support rods | 39 |
| 3.2 Mounting of the instrument to a stand | 39 |
| 3.3 Mounting the stirring paddles | 40 |
| 3.4 Mounting the shaft protector (optional) | 40 |
| 4. Setup and operation | 41-44 |
| 4.1 Intended use | 41 |
| 4.2 The front panel | 41 |
| 4.3 Description of the front panel functions | 42 |
| 4.4 Description of the rear panel | 43 |
| 4.5 How to display the system info | 44 |
| 5. Safety functions | 44-45 |
| 5.1 Rotor stuck protection | 44 |
| 5.2 Overload protection | 44 |
| 5.3 Overheating protection | 44 |
| 5.4 I/O connector (9-pin SUB-D) | 45 |
| 6. The RS232 interface | 45-51 |
| 6.1 Description of the RS232 interface | 45, 46 |
| 6.2 RS232 connection between PC and stirrer | 46, 47 |
| 6.3 Format of an RS232 command | 47 |
| 6.4 Format of the controller handshake | 47, 48 |
| 6.5 RS232 commands of the stirrer | 48, 49 |
| 6.6 Programming examples | 49 |
| 6.6.1 Writing parameters to the stirrer | 49, 50 |
| 6.6.2 Reading parameters from the stirrer | 50 |
| 6.7 Connecting more than one unit to a computer | 51 |
| 7. The setup menu | 52 |
| 7.1 Description of the setup menu | 52, 53 |
| 7.2 Example: how to change setup settings | 53, 54, 55 |
| 8. Error messages | 55 |
| 9. Cleaning and maintenance | 56 |

| | |
|---|--------|
| 10. Dismantling, transport and storage | 56-57 |
| 10.1 Dismantling | 56 |
| 10.2 Transport and storage | 56, 57 |
| 11. Technical data | 57, 58 |
| 12. Repairs | 58 |
| 13. Warranty | 59 |
| 14. Liability | 59 |
| 15. Disposal | 59 |

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Laboratory Stirrer R100 C

1. Suitable for stirring media with medium to high viscosity
2. Display of the set and real speed
3. User programmable limits for speed and torque
4. Motor speed 40-2000 rpm
5. max. torque 150 Ncm
6. Timer function
7. Hollow shaft for 1-10,5 mm diameter paddles
8. Overheating protection
9. RS232 interface

For the safe use of the stirrer we recommend the use of the following accessories:
U-stand and special clamp

1. General information

1.1 Unpacking the instrument

The stirrer R100 C was manufactured and tested according to DIN EN 61010. According to these regulations, the units are designed to meet the requirements for safe and correct operations. To ensure the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred during transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

1.2 Scope of delivery

Your delivery should contain the following:

- 1 stirrer R100 C
- 2 support rods
- 1 chuck key
- 1 instruction manual

Stirrers and a stand are not included in the scope of delivery. For the safe use of the stirrer we recommend the use of the following accessories: U-stand and special clamp (see 1.3)

Note:

If the instrument is not damaged and all parts are complete, you may start to operate the device after reading the instruction manual.

Important note:

Read this instruction manual carefully before operating the instrument. Please store the manual in a place easily accessible for every user. Should there be any further questions after reading the manual concerning the setup, operation or warranty, please contact:

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1.3 Accessories (optional)

| | |
|--|------------------------|
| U-stand (recommended) | Marienfeld No. 7205225 |
| Special clamp for mounting the instrument to the stand (recommended) | Marienfeld No. 7205226 |
| Stirring shaft protector | Marienfeld No. 7205227 |
| Jumbo-clamp to mount the shaft protector to the stand | Marienfeld No. 7205228 |
| RS 232 interface cable D-Sub 9-pin 1,8 m | Marienfeld No. 7205229 |
| USB to RS232 adapter | Marienfeld No. 7205230 |
| Software Stirrer Control | on request |

1.4 Important instructions

A SCHUKO plug (DIN 49441, 16 A, 250 V) is normally shipped with our instruments. North America instruments are shipped with standard US plugs (NEMA Pub. No. WDI1961 ASA C 73.1. 1961 page 8 15A 125V).

When operating the instruments in countries with different AC plug systems, use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation. Ensure that your local supply voltage matches the indication on the instrument's specification plate!

The instrument is earthed. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

2. Safety instructions

2.1 Description of safety symbols in the instruction manual

This symbol distinguishes notices which have to be strictly observed by the user. Disregard of warnings may result in the impairment of serviceability as well as in physical harm to the user.

This symbol distinguishes interdictions which have to be strictly observed by the user. Disregard of warnings may result in the impairment of serviceability as well as in physical harm to the user.

This symbol distinguishes notices which have to be strictly observed by the user to ensure safe operation of the unit.

Caution
Danger of fire or explosion!

Note regarding repair or maintenance

Note regarding the mains cable connection

2.2 Description of safety symbols on the instrument

Warning against general danger:

This symbol indicates that it is imperative to read and understand the instruction manual prior to operating the instrument. Please highlight points which require special attention in your field of application so they are not overlooked. Disregard of warnings may result in impairment of serviceability as well as in physical harm to the user.

2.3 General safety instructions

Caution:

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

Caution

Please do not use this instrument in insecure environment especially not in explosive ambience. Danger of life!

**Caution:**

Instructed users only may operate the instrument.

**Caution:**

When connecting the instrument to an AC power outlet, ensure that your local supply voltage corresponds to the specifications indicated on the instrument.

**Caution:**

Please act with caution when working in the vicinity of flammable or explosive substances and observe safety data sheets. This instrument is not explosion-proof.

The On/Off switch does not disconnect the instrument from the power source. Remove the plug from the AC power outlet to disconnect the instrument from the mains supply entirely.

**Note:**

Make sure that all safety information on the instrument is clearly visible during operation.

**Caution:**

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

Connect the instrument to an earthed AC power outlet only.

**Caution:**


Make sure the support used is strong enough to hold the instrument and is standing on a solid and horizontal surface.

3. Mounting

3.1 Mounting of the support rods

chuck key (B)


support rods (A)

Screw the two support rods (A) into the threaded holes at the back of the overhead stirrer (see chapter 4.4). Use the chuck key (B) to firmly tighten them.

3.2 Mounting the instrument to the stand

special clamp (C)

knurled screw (E)


support rod (D)

Mount the special clamp (C) to the rod of the stand (D). Insert the support rods (A) of the stirrer into the openings of the special clamp (C) and turn the knurled screws (E) until the drive unit is firmly attached to the stand.

The position of the stirrer can be changed once the motor is mounted. Hold the stirrer firmly, loosen the screws and bring the stirrer into the desired position. Hold the stirrer with one hand and fasten the screws again with the other hand.


3.3 Mounting the stirring paddle

Insert the stirring paddle (G) as shown in the picture into the chuck (F) and firmly tighten it in the desired position using the chuck key (B). The immersion depth of the paddle can be adjusted by loosening the chuck and axially moving the stirring tool.

3.4 Mounting the shaft protector (optional accessory)

We recommend to use the shaft protector for secure working with the stirrer (optional accessory Marienfeld no. 7205227).

Attach the fastening rods (H) to the stand rod in front of the shaft protector using for example the Jumbo clamp (Marienfeld no. 7205228).

Position the bores (K) of the protector next to the screws (J), insert the screws (J) into the bores and turn the protector to the right as far as it will go.

The screw (L) is used to adjust the overall height of the protector.

Check that the stirring shaft protector is held in position securely prior to each use and also at regular intervals.

The position of the protector must be adjusted only when the stirrer is turned off and the power supply is disconnected.

4. Setup and operation


4.1 Intended use

The stirrer R100 C is designed for the skilled use in laboratories. It is suitable for stirring and mixing of liquids with medium to high viscosity. For correct use the instrument must be fixed to an appropriate stand.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and make sure that the stirrer is not exposed to a corrosive atmosphere.

4.2 The front panel

The innovative key-pad facilitates data entry and incorporates a two-line alpha-numeric LCD-display (1) which shows all relevant data. The entry and changing of timer or setup settings is made via the setting control knob (5), the motor speed is set via the right speed control knob (7).


4.3 Description of the front panel functions

| | |
|----------------------------------|---|
| LCD-display (1) | The LCD display shows all relevant system data |
| Timer LED (2) | indicates that the timer is active and the unit will be automatically switched off when the time has expired |
| Timer button (3) | If required, pressing this key initiates the timer function to shut down the stirrer. To change the timer setting, use the setting control knob. The time will be shown as days : hours : minutes |
| Select Display button (4) | Press this key to step forward from one display option to the next one. Press and hold the key for at least 3 seconds to zero the relative torque measurement (Mr). Press and hold the key for at least 5 seconds to switch back from the relative torque measurement (Mr) to the absolute torque measurement (Ma). The same applies to the power output measurement function. |
| Setting control knob (5) | entry and change of timer or setup settings |
| Load LED (6) | indicates that the unit is not able to hold the set motor speed due to a too high torque of the motor shaft. The instrument now will reduce the current motor speed so that the maximum allowed motor torque is not exceeded. |
| Speed control knob (7) | setting the speed manually |
| Power LED (8) | indicates that the unit is switched on |
| Power on/off key (9) | Press this key to turn off the instrument. |


4.4 Description of the rear panel

On the rear panel of the R100 C you will find 2 threads (a) to mount the support rods (A), a 9-pin SUB-D connector (10), a fan (11) and the mains power supply cable (12)

Support rod (A)

| | |
|--|--|
| <p>9-pin SUB-D connector (10)</p> | <p>Input/output for:</p> <ul style="list-style-type: none"> • RS232 interface (full duplex) • analogue input for speed setting via voltage signal input range: 0-10 V, operation range programmable • analogue output for speed, torque or power monitoring output range: 0-4 V, operation range programmable |
| <p>Mains power supply (12)</p> | <p>Make sure the voltage printed on the R100 C corresponds to the voltage of your mains!</p> |

4.5 How to display the system info

For additional user information the instrument tracks the total operation time and the frequency of switching on and off.

Press and hold the timer key (3) whilst switching on the stirrer. The instrument initialises and displays the system info screen. Press any key to display the next information.

The system info is shown on the LCD-display (1) and comprises the following details:

1. total operation time of the unit (days : hours : minutes)
2. frequency of switching on and off

5. Safety functions

5.1 Rotor stuck protection

In case that the torque exceeds a certain limit over the time, the instrument will assume that the stirrer (rotor) was suddenly blocked and the instrument is switched off immediately. The cut-off limit sensitivity can be adjusted in the setup menu (→ Cut off limit)

5.2 Overload protection

If the torque or power drawn from the stirring motor exceeds the allowed range, the unit will automatically reduce the motor speed until the torque is reduced to the admissible value. The red LED load (6) indicates this incident.

As long as the Load LED (6) lights up, the motor can't hold the set motor speed.

5.3 Overheating protection

In case that the temperature of the motor exceeds the maximum operating temperature, the unit will automatically reduce the power output (speed and torque) in order to reduce the temperature to a safe operation point (= fold back characteristic). This protects the stirrer against damage by constant overload operation. The overload operation is also displayed by the red Load LED (6).

A cooling fan (11) on the rear side of the stirrer is automatically switched on additionally.

As soon as the temperature of the motor exceeds 70 °C the instrument is switched off automatically. The safety feature is fully redundant.

Important note:

Before re-starting the unit after cooling down, switch off at first and after a while switch it on again using the on/off button.


5.4 I/O connector (9-pin SUB-D)

Pinout of the I/O connector on the rear side of the instrument:

- RS232 control signals
- Analogue input (0-10 V) for speed override
- Analogue output (0-4 V) for recording of process parameters (speed, torque or power)

Pinout of the I/O connector:

| Pin number | Description |
|------------|-------------------------|
| 1 | +5 V (max. 10 mA) |
| 2 | TxD (transmit data) |
| 3 | RxD (receive data) |
| 4 | n.c. |
| 5 | GND (signal ground) |
| 6 | analogue output (0-4 V) |
| 7 | n.c. |
| 8 | n.c. |
| 9 | analogue input (0-10 V) |

6. The RS232 interface

6.1 Description of the RS232 interface

The RS232 interface enables remote control of all stirrer functions. In this configuration a computer can access each connected unit via a sub address (daisy chaining). This address can be programmable for every unit (→ setup menu, slave-address).

For daisy-chaining the transmitting line (TxD) of the PC is connected to the receiving line (RxD) of the first unit. The transmitting line of this unit then is connected to the receiving line of the next unit in the chain. The transmitting line of the last unit in the chain then is connected to the receiving line of the computer which closes the chain (ring).

For addressing a specific unit in a daisy-chain, each controller carries a so called slave-address which can be any number from 1 to 255. The default slave address is 1. The slave-address 0 is defined as a general call-address to which all controllers will respond. The slave-address can be set or changed in the setup menu.


All instruments in a chain have to be adjusted to the same baud rate.

A software to control a connected unit is available on request.

Pinout of the RS232 connector:

| Description | pin number |
|---------------------|------------|
| TxD (transmit data) | 2 |
| RxD (receive data) | 3 |
| GND (GROUND) | 5 |


Parameters of the serial interface (1200-4800,8,N,1):

Baud rate: 1200-4800 Baud (1200 default, see „Setup menu“
 Databits: 8 Bit
 Parity: none
 Stopbits: 1


6.2 RS232 connection between PC and stirrer

a) PC with 9-pin RS232 plug

b) PC with 25-pin RS232 plug

6.3 Format of an RS232 command

Each command sent to the instrument must have the following format:

ADR, CMDCODE, PARAMETERLIST <CR>

Description:

ADR: slave address of the unit to execute the command
 CMDCODE: command code
 PARAMETERLIST: 1 to 6 parameters separated by comma
 CR: the command string must be terminated by Carriage/Return (ASCII 13)

6.4 Format of the controller handshake

After receiving an RSR-command the unit will:

3. send the received command to the next unit (or back to the PC – daisy chaining)
4. answer with a handshake string which is defined as follows:
 ADR, „HS“, RETCODE, PARAMETERLIST CR

Explanation:

ADR: Slave address of the unit sending the handshake
 RETCODE: error code (see table on next page)
 PARAMETERLIST: 1 to 6 parameters separated by comma (see table 2, command codes)
 CR: the handshake as any command is terminated by ASCII-code 13 (CR)

Table 1

| Return Code | Explanation | Parameterlist |
|-------------|---|-----------------------|
| OK | command executed, no error | see table 1 |
| UC | unknown command | none |
| PA | wrong parameter number (too few or too many parameters specified) | none |
| NA | command is not allowed in actual operation mode | actual operation mode |
| PR | at least one parameter is out of range | none |
| PL | at least one parameter is too long | none |
| DF | unknown data format | none |

6.5 RS232 commands of the R100 C stirrer

| CMD. Code | Explanation | Parameter list | Range |
|-----------|--|--|-----------------|
| RAC | Read actual settings of the unit | 1. Dummy parameter to initiate transfer → controller sends handshake A) actual motor speed (rpm) B) actual torque (relative) C) actual torque (absolute) D) actual power (relative) E) actual power (absolute) | 1 |
| WTR | Write timer | A) Timer value in [sec] | 0..2592000 |
| RTR | Read timer | 1. Dummy parameter to initiate transfer → controller sends handshake: A) timer value (seconds) → 0 = Timer disabled | 1 0..2592000 |
| WSE | Write motor (SET) define motor speed setpoint in rpm | A) Motor speed in rpm 0 → motor switched off (speed =0) If selected speed exceeds the permissible limit, the value will be reduced to the acceptable limit (→ Setup menu, max. MotorRpm) | 0..2000 |
| RSE | Read motor setpoint (SET) | 1. Dummy parameter to initiate transfer → controller sends handshake: A= motor setpoint in rpm | 1 0..2000 |
| RSS | Read system status | 1. Dummy parameter to initiate transfer → controller sends handshake: 1. actual system status code 1 → normal operation 2. overload condition 0 → no overload condition 1 → motor is under overload operation | 1 1 0/1 |


| | | | |
|-----|--|---|---|
| RTY | Read type information of unit | 1. Dummy parameter to initiate transfer → controller sends handshake: A) name of the connected device B) firmware version number C) ON/OFF counts D) total operation period in minutes | 1 string string number number |
| WSM | Keyboard or RS232 speed control (write serial mode). Select if setting for motor speed is accepted from front panel (speed control knob) or via RS232 command (WSE) | A) Control parameter 0 → front panel control of motor speed 1 → RS232 control of motor speed (→ via WSE command) | 0/1 |
| KEN | Lock/unlock front panel keyboard | 1. lock/unlock 0 → front panel keyboard is disabled (except On/Off button) 1 → front panel keyboard is enabled | 0/1 |
| WSA | Set RS232 slave-address + renumber slaves | 1. New slave-address of unit | 1...255 |
| OFF | Switch off unit | 1 Safety parameter | 1234 |
| WEE | Save setup settings to EEPROM | 1. Safety parameter | 4321 |

6.6 Programming examples

The following application examples show how RS232 commands can be sent to the stirrer and how actual parameters can be read out from the stirrer.

6.6.1 Writing parameters to the stirrer

If, for example, the motor speed should be set to 500 rpm the following two command strings have to be sent to the unit with slave-address 1:

1,WSM,1 ; 1. command: disable front panel control and enable RS232 control
1, WMS,500 ; 2. handshake (response) to the 1. command (= accepted, no error)

The R100 C controller now sends the following strings to the next unit (or back to the PC if no further instruments are connected):

1,WSM,1 ; echo of the 1. command
1,HS,OK ; handshake (response) to the 1. command (= accepted, no error)

1, WMS,500 ; echo of the 2. command
; handshake (response) to the 2. command (= accepted, no error)

6.6.2 Reading parameters from the stirrer

To read the actual values of motor speed, torque and power output send the following command to the unit with slave-address 1

slave-address = 1
command code (RAC = read actual settings)
1. parameter
here: only dummy parameter to fulfil protocol
(at least one parameter must be sent)

1,RAC,1 ; see table 2: RAC-command

The controller now sends the following two strings to the next unit (or back to PC):

1,RAC,1 ; echo of the received command
1,HS,OK, 300.0, 50.0, 56.0, 15.7, 17.59 ; handshake (response) with return parameters

absolute power output (Pa) in Watt
relative power output (Pr) in Watt
absolute torque (Ma) in Ncm
relative torque (Mr) in Ncm
1. parameter is actual motor speed in rpm
return code (see table 1), OK = no error
describes that this is a handshake and not a command (handshake designator)
slave-address of the unit which sent the handshake

Note:

every command or handshake sequence is terminated by ASCII-Code 13 (CR)


6.7 Connecting more than one unit to a computer

The powerful protocol of the RS232 interface allows also for connection of other laboratory instruments like hotplate stirrers, dispensers etc. to one RS232 interface only on the computer side. By assigning different slave-addresses to each connected unit, the computer can access each unit individually.

For connection of more than one unit to a computer a special Y-cable is used for every additional unit (RS232-Y-cable, available on request).

For application adaptation a Windows® program is available on request.

7. The setup menu

7.1 Description of the setup menu

Global settings are stored in a non-volatile memory (EEPROM). To modify these global settings, press and hold the "Select Display Key" (4) while switching the unit on with the Power On/Off key (9).

The following settings can be adjusted in the setup menu of the R100 C:

| | |
|----------------|--|
| RS232 address | Defines the address under which the stirrer can be accessed using RS232 commands. factory setting: 1 range: 1-255 |
| Baud rate | Baud rate of the RS232 interface factory setting: 1200 Baud options: 1200, 2400, 4800 Baud |
| Max MotorRpm | User selectable maximum motor speed in rpm. factory setting: 2000 rpm range: 40-2000 rpm |
| Max torque | max. permissible torque on shaft in Ncm factory setting: 150 Ncm range: 10-150 Ncm |
| CutOff limit | Cut-off limit sensitivity in percent (→ rotor stuck detection) If the torque increase over the time exceeds a certain limit, the unit will assume that the rotor was suddenly blocked. If this happens the unit will shut down immediately. factory setting: 60% range: 10-100% (100% disables this safety check) Note: The lower the entered value, the more sensitively the unit will react on sudden torque increases. |
| Analogue input | Allow or block speed setting via external analogue voltage control factory setting: disabled options: enable/disable |
| An Inp. F=0 | Defines the control voltage (→ analogue input), at which the motor has to be stopped (speed = 0). The voltage at the analogue input functions as multiplication factor onto the frontally set speed (→ speed control knob). This value reflects the input voltage at which the multiplication factor is zero. This setting is only relevant if the analogue input is enabled (see above). factory setting: 0.0 V range: -20.0 V to +30.0 V (<u>only to define the load line, max. applied voltage: 0-10 V</u>) |

| | |
|--------------|---|
| An Inp. F=1 | This value reflects the input voltage at which the multiplication factor of the speed is 1. This setting is relevant only if the analogue input is enabled. factory setting: 5.0 V (50 % of the torque set on the instrument; 10 V is 100 % of the torque set) range: -20.0 V to + 30.0 V (<u>only to define the load line, max. applied voltage: 0-10 V</u>) |
| An Output | Selects which function/value has to be monitored on the analogue output pin. factory setting: DISABLED Available options: Torque-R (relative torque) Torque-A (absolute torque) Power-R (relative power) Power-A (absolute power) DISABLED (output de-activated, 0V output) |
| An. Out U=0V | This value defines the operation value (speed, torque or power) which is related to a voltage outlet of 0.0V at the analogue pin output (pin no. 6) |
| An. Out U=4V | This value indicates which operation value (speed, torque or power) is related to the output voltage of 4.0V at the analogue pin output (pin no. 6). |

7.2 Example: How to change setup settings

Example for the necessary steps to set or define the maximum permissible motor speed in the setup menu:


- Press and hold the "Select Display" button (4) whilst switching on the instrument. The instrument initialises and switches into the setup menu

Note:

The „Select Display“ key (4) must be pressed down until the setup menu appears on the display (1)

- Turn the handwheel (5) until „Max MotorRpm“ is shown on the display (1)
- Press the Timer button (3), to change the max. motor speed (→ edit arrow is shown)
- Turn the handwheel (5), to adjust the value
- Press the „Select Display“ key (4), to quit the setup menu
- Press once again the „Select Display“ key (4) to save the setup parameters

**Note:**

The different options/parameters in the setup menu can be displayed/selected by simply turning the handwheel (5). To change a displayed value, press the Timer key (3) and an edit arrow (→) will appear in front of the parameter to be changed. As long as the edit arrow is shown on the display, turning the handwheel (5) will change the displayed value. If the handwheel (5) is not turned for approx. 3 seconds, the edit arrow will disappear and another menu option can be selected by turning the handwheel (5).

When editing the setup parameters, the "Select Display" button (4) is used to leave the setup menu, or to confirm a selection. The "Timer" button (3) is used to edit/change a displayed value.

The setup menu can only be left by pressing the Power On/Off key (9) or the "Select Display" key (4).

Pressing the Power On/Off key (9) in the setup menu will prevent the unit to store eventually changed settings to non-volatile memory. After pressing the "Select Display" key (4) the unit will ask if the settings should be stored to non-volatile memory. To confirm this, press "Select Display" (4) again. If the parameters should not be stored, turn the handwheel (5) until "NO" is displayed and then press the "Select Display" button (4).

8. Error messages

The R100 C stirrer disposes of an on-board self-diagnostic program that detects possible failures and in case of such a failure, the instrument shuts down. The reason for this shut down is displayed for about 3 seconds when switching on the stirrer again.

| Error message: | Reason for shut down: |
|----------------|--|
| Time-out | The time set has expired (timer key) |
| rotor stuck | A too fast increase of the rotor torque was detected (see: setup menu) |
| Overtemp. | The motor temperature was too high |
| RS232 Off | The unit was switched off by an RS232 command |
| Sensor Err. | The temperature measurement sensor or circuit is damaged |

9. Cleaning and maintenance

The outer casing consists of painted or anodised aluminium and a chemically resistant splash-proof membrane key pad. To clean the instrument, use a damp cloth and a mild soap solution.

Do not use steel wool or any aggressive cleaning solution to clean the instrument. This could cause damages of the casing.

Do not attempt to repair the instrument in case of mal function. There are no user-serviceable parts in the instrument. **Any attempt by the user to repair the unit will immediately render the warranty null and void!**

10. Dismantling, transport and storage

10.1 Dismantling

Switch off the instrument and disconnect from the mains. Remove any glass beakers and any other equipment around the stirrer.

chuck (F)
chuck key (B)
stirring paddle (G)

Remove the stirring paddle (G) by loosening the chuck (F) with the chuck key (B). Loosen the special clamp (C) and remove the instrument. Unscrew the support rods (A) (see chapter 3.2)

10.2 Transport and storage

Place the instrument in its original packaging if possible. Each module should be protected by air bubble film against shocks from outside. When placing the individual modules in the carton, make sure there is adequate interspace between the parts which should be filled with padding material.

In case you do not use the original packaging, please mark your package with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)

For an appropriate storage please observe the following:

4. Dismantle all parts from the instrument
5. Each part should be packed separately in air-tight plastic bags
6. Observe the storage ambient:
 - max. ambient temperature: room temperature up to + 40 °C
 - max. humidity: 80 %

11. Technical data

| | |
|------------------------------------|---|
| Power | 115/230 V, 300 W, 50/60 Hz (see rear panel) |
| LCD-Display | bright back-lit LCD display, showing all relevant system data. Simultaneous display of set and real values of motor speed, torque, or power output. |
| Handwheel/setting control knob (5) | allows fast and intuitive adjustment of timer and system settings |
| Timer function | programmable switch-off timer, programmable from 1 minute to 99 days. Resolution: 1 minute |
| Torque output display | displays true torque at shaft in Ncm (relative or absolute readout possible) Accuracy: 5% |
| Power output display | displays true power output at shaft in watts. (relative or absolute readout possible) Accuracy: 5% |
| Motor specification | <ul style="list-style-type: none"> • torque microprocessor controlled for easy operation and accuracy • soft-start and soft-stop logic prevents splashing of liquid • precise adjustable motor speed from 50 to 2000 rpm • in steps of 10 rpm • max. 220 watts shaft power output • extra long life, high performance, brush motor • silent belt drive power transmission • max. torque: 150 Ncm (20 min) • continuous torque output (at 1100 rpm): 100 Ncm • maximum limitation of torque (programmable) |
| Multistep safety system | detects and protects against the following dangerous situations: <ul style="list-style-type: none"> • rotor stuck/block detection • motor overheating monitoring • selftest of all safety functions after switching on |
| Auxiliary functions | User programmable limits for: <ul style="list-style-type: none"> • maximum motor speed (40-2000 rpm) • maximum allowed torque at shaft (10-150 Ncm) |
| RS232 interface | programmable Baud rate: 1200,2400,4800; 8,N,1 Full remote access of all stirrer functions (read and write) daisy-chain feature by sub addressing allows to connect up to 255 units to only one serial interface on the computer side. Windows® software optional available |
| Analogue input to write: | Range: 0-10 V Programmable range Allows speed setting via external analogue voltage control |

| | |
|-----------------------------------|--|
| Analogue output to read: | Range: 0-4 V Range of output and control (speed, torque and output) is programmable |
| System info | For additional user information total operation time and total frequency of switch on/off can be retrieved (see chapter 4.5) |
| Permissible ambient temperature: | 5-40 °C within operating area, e.g. fume cupboard |
| Max. permissible humidity: | 80 % RH |
| Safety class acc. to DIN EN 60529 | IP42 |
| Case dimensions: | 78x200x230 mm (W x D x H) |
| Weight: | 4.7 kg |

12. Repairs

Attention:

Do not attempt to repair the instrument in case of mal function. There are no user-serviceable parts in the instrument. **Any attempt by the user to repair the unit will immediately render the warranty null and void!**

When returning instruments for repair which have come into contact with hazardous substances, please:

- provide precise information on the relevant medium
- take protective measures to ensure the safety for our receiving and maintenance personnel
- mark the package as appropriate for hazardous materials

13. Warranty

Paul Marienfeld GmbH & Co. KG provides a guarantee for a period of 24 months from date of delivery for the instrument purchased by you. This warranty applies only to defects in material or workmanship. There is no warranty for defects or malfunctions caused by neglect, improper use or improper service.

In case of defects in material or workmanship, the instrument will be repaired or replaced. In case of replacement, the warranty for the new unit will persist for the remaining period of the primary warranty period. Charges for returning the instrument have to be borne by the customer.

Any claims regarding guarantee and liability are limited to professional and proper use of the instrument and observance of the general safety instructions and this instruction manual.

14. Liability

Paul Marienfeld GmbH & Co. KG is not liable for any material or personal damages caused by improper use and does not accept liability for consequential damages.

It is under the responsibility of the user to find out whether the device is suitable for his application. If in doubt, clarify this with your dealer or with Paul Marienfeld GmbH & Co. KG directly.

The warranty expires in case of removal or alteration of the identification label and/or the serial number of the instrument.

15. Disposal

Damaged and/or waste electric or electronic equipment has to be disposed of at the designated Recycling Depot.

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